

# Pocket Spot Checks at All Times

## iM3s Vital Signs Monitor



iM3s  
Vital Signs Monitor



# Small Outside, Smart Inside

Being the smallest vital signs monitor among its kind, iM3s presents a fashionable handheld design, delivering extraordinary portability throughout spot checks. Though small in size, it integrates all necessary vital measurements together, including SpO<sub>2</sub>, NIBP and optional TEMP.

## A Dedicated Spot Monitor

With its pleasing size and large capacitive touch screen, iM3s could be operated by one hand, bringing comfortable and familiar user experience to daily applications. For each test, all measurements and custom parameters could be displayed on the same screen, while the test results being saved automatically in its built-in non-volatile memory or transfer to hospital information system.



300g  
Lightweight Design



IP44  
Waterproof & Dustproof



Shock Resistance



≥6 hours  
Minimum Run Time



Wi-Fi  
Network Connection



e-Link  
Wireless Connection



20,000 sets  
Internal Storage



Barcode Scanner




## Tabletop Expansion

Plugged onto an optional tabletop dock, the iM3s easily fits into regular tabletop applications, enabling more possibilities, including AC charging, optional thermal recording, optional temporal temperature measurements, etc.



Hand Carry Design



AC Charging



Option  
Temporal TEMP



Option  
Thermal Recorder

## Powerful functions with flexible configuration

### Main Unit

- Standard Configuration: SpO<sub>2</sub>, NIBP, Capacitive touch screen
- Optional Configuration: Li-ion battery, Built-in barcode scanner, 2.4G/5G Hz Wi-Fi, e-Link

### CS-04 Charger Stand

- AC power supply

### CS-05 Charger Stand

- Standard Configuration: Built-in Li-ion battery, AC power supply, USB port, RJ45 port
- Optional Configuration: Thermal recorder, Exergen TEMP (TAT-5000S), Quick TEMP(F3000)